

Ultralow temperature and Cryogenic freezer

The representative product of OPERON Co., Ltd, -156°C cryogenic freezer for preserving stem cell is the compression freezer which is used to preserve samples for a long time without deterioration and does never use liquid nitrogen. OPERON cryogenic freezer which overcame high cost, high risk, inconvenience and sample contamination occurring when using liquid nitrogen is registered to the international patent (the United States, Germany, China and Korea).

OPERON cryogenic freezer is applied with cryogenic cooling system, and it is an innovative cryogenic product with electric method replacing liquid nitrogen tank. The temperature range in which ice is not formed is called "AMORPHOUS", and the approximate range is $-135^{\circ}\text{C} \sim -150^{\circ}\text{C}$. If the samples go through this temperature range rapidly, it is able to preserve the samples without deterioration from several years to maximum several decades because ice crystals are not formed. Other than this, we have various products like independent control type freezer with separate chambers, shelf-cooling type freezer which has powerful freezing and temperature uniformity, power plus freezer which is appropriate for surrounding temperature $+40^{\circ}\text{C}$ and ultra safe (Twin Heart) freezer which has 2 compressors which maintain -90°C during the operation, and when one compressor is broken down, another compressor maintains -80°C .



- Cryogenic freezer (CFQ)
- Ultralow-upright (DFU)
- Chest (DFC)
- Double door (DFUD)
- Shelf-cooling type (DFU-V/DFUD-V)
- Power plus (DFU-P/DFUD-P)
- Ultra safe-Twin Heart(DFU-CUS/DFC-CUS)

Cryogenic Freezer-CFQ


Cryo Freezer conqueror -156°C

Product features and specifications

- The world first freezer which makes maximum cryogenic -156°C without liquid nitrogen
- Application of cryogenic cooling system (Auto cascade system) registered to the international patent
- Use of CFC-free special mixed refrigerants which OPERON Co., Ltd developed
- Application of aerodynamically designed pan wings equipped with huge capacity and high performance
- Hermetic compressor for LBP with excellent performance
- CE certification product for preserving cell
- Non scaling stainless chamber and round corner design for easy cleaning
- Use of 200mm high density urethane insulating materials to minimize cold air loss
- Chamber mold and double door with elegant design
- Strong structure frame and impact resistant powder coating finishing
- Filter cover for easy cleaning
- Chest type cabinet structure
- Temperature for hard condition resistant test +32°C
- Temperature measuring sensor: platinum PT-100Ω(Class A 0.15 grade)
- Prevention of frost with Rim Heating system
- RS-232 Data port and Interface which is available for store, control and monitor temperature data via user's PC






| Items | OPERON ELECTRIC FREEZER | LIQUID NITROGEN FREEZER |
|---------------------|---|--|
| Economic efficiency | It is available only with the lower cost electronic connection. | Expensive liquid nitrogen tank and the liquid nitrogen should be supplemented consistently. |
| Safety | It is safer from pollutants. When blackout, it secures safety with BACK UP SYSTEM if the temperature goes up. | It has a risk to be contaminated by bacteria or fungi inside the tank |
| Accuracy | Temperature difference between the upper part and the lower part inside is $\pm 5^{\circ}\text{C}$, so the temperature is more even. | The temperature of the lower part inside the tank is -196°C , and the upper part is -70°C . The temperature difference is huge. This has a significant impact on viability of stem cell. |
| Convenience | Internal space can be used efficiently for keeping samples, and it is simple to perform work. | Rack for keeping samples is very inconvenient for sample work because of the liquid nitrogen tank's structure, and it is difficult to use the space inside the tank efficiently. |



Stem cell needs long-term preservation for 10~20 years.

It's been stated on various papers that cell preservation using liquid nitrogen has lots of problems. OPERON products resolved the above problems perfectly.



[References for reliability evaluation of cryogenic technology research center affiliated to OPERON]

Control system

- Micro processor control system /Digital temperature indicator which adjust the temperature by 0.1°C
- Functions for setting double passwords to prevent access or control except administrators
- Data locking function

Alarm system and safety devices

- Alarm for filter cleaning period (Warm Condenser alarm)
- Alarm for high/low temperature of sight and auditory
- Equipped with pressure switch
- Equipped with overheating prevention switch
- Auto return system after blackout

Cryogenic Freezer (Cryo freezer)

| Model | CFQ-150 | CFQ-152 | CFQ-156 | CFQ-232 | CFQ-300 |
|--------------------|--|---|--|---|---|
| External size(mm) | W1614 x D890 x H1105 | W1614 x D890 x H1105 | W1614 x D890 x H1105 | W1676 x D890 x H1105 | W2170 x D900 x H1065 |
| Internal size(mm) | W600 x D400 x H600 | W600 x D400 x H600 | W600 x D400 x H600 | W770 x D460 x H665 | W1050 x D470 x H610 |
| Range for usage | $-90^{\circ}\text{C} \sim -150^{\circ}\text{C}$ | $-95^{\circ}\text{C} \sim -152^{\circ}\text{C}$ | $-100^{\circ}\text{C} \sim -156^{\circ}\text{C}$ | $-90^{\circ}\text{C} \sim -156^{\circ}\text{C}$ | $-80^{\circ}\text{C} \sim -140^{\circ}\text{C}$ |
| Inside measurement | 144L | | | 235L | 300L |
| Electric capacity | Capacity necessary for installment AC220V 380V 3ph (50/60Hz) | | | | |
| Weight | 380Kg | | | | 450Kg |

Options

Rack / Box / LN2 Backup system / 7"LCD touch screen and USB (3 years Data) - Basic -1channel (sensor 3channel, 6channel options) / U-system (remote alarm system and wireless control system – telecommunication subscription needed / Automatic voltage regulator / Cryogenic gloves / Battery backup (Control panel and alarm) / Recorder / VIP

Ultralow temperature freezer-DFC

Chest Type

Product features and specifications

- Realization of powerful and stable -86°C freezing performance with application of Two stage Auto cascade system
- Use of CFC-free special mixed refrigerants
- Application of aerodynamically designed pan wings equipped with huge capacity and high performance
- Hermetic compressor for LBP with excellent performance
- Use of 130mm high density urethane insulating materials to minimize cold air loss
- Non scaling stainless chamber
- Strong structure frame and impact resistant powder coating finishing
- Round corner design for easy cleaning
- Prevention of condensation and the chill loss with two-layered packing
- Filter cover for easy cleaning and washable filter
- Chest type cabinet structure
- Temperature for hard condition resistant test +32°C
- Temperature measuring sensor: platinum PT-100Ω(Class A 0.15 grade)
- Prevention of frost with Rim Heating system
- RS-232 Data port and Interface which is available for store, control and monitor temperature data via user's PC



Control system

- Micro processor control system /Digital temperature indicator which adjust the temperature by 0.1°C
- Functions for setting double passwords to prevent access or control except administrators
- Data locking function

Alarm system and safety devices

- Alarm for filter cleaning period (Warm Condenser alarm)
- Alarm for high/low temperature of sight and auditory
- Equipped with pressure switch
- Equipped with overheating prevention switch
- Auto return system after blackout

Ultralow temperature freezer (Chest freezer)

| Model | DFC-37 | | | DFC-84 | | | DFC-200 | | | DFC-300 | | | DFC-400 | | | DFC-500 | | | DFC-600 | | |
|--------------------|---|-------|-------|---------------------|-------|-------|----------------------|-------|-------|----------------------|-------|-------|----------------------|-------|-------|----------------------|-------|-------|----------------------|-------|-------|
| Temperature range | AE | BE | CE | AE | BE | CE | AE | BE | CE | AE | BE | CE | AE | BE | CE | AE | BE | CE | AE | BE | CE |
| | -40°C | -55°C | -86°C | -40°C | -55°C | -86°C | -40°C | -55°C | -86°C | -40°C | -55°C | -86°C | -40°C | -55°C | -86°C | -40°C | -55°C | -86°C | -40°C | -55°C | -86°C |
| External size | W430 x D560 x H920 | | | W700 x D653 x H1125 | | | W1365 x D720 x H1101 | | | W1665 x D720 x H1101 | | | W1965 x D720 x H1101 | | | W2265 x D720 x H1101 | | | W2565 x D720 x H1101 | | |
| Internal size | W235 x D365 x H425 | | | W437 x D387 x H500 | | | W600 x D470 x H710 | | | W900 x D470 x H710 | | | W1200 x D470 x H710 | | | W1500 x D470 x H710 | | | W1800 x D470 x H710 | | |
| Inside measurement | 37L | | | 84L | | | 200L | | | 300L | | | 400L | | | 500L | | | 600L | | |
| Electric capacity | Capacity necessary for installment AC220V 1ph (50/60Hz) | | | | | | | | | | | | | | | | | | | | |
| Weight | 40kg | | | 150Kg | | | 255Kg | | | 290Kg | | | 315Kg | | | 340Kg | | | 365Kg | | |

Ultralow temperature freezer (Power Plus model-Chest type)

- Chest type Power Plus model- the freezer appropriate for the surrounding temperature +40°C and over

| Model | DFC-200P | | | DFC-300P | | | DFC-400P | | | DFC-500P | | | DFC-600P | | |
|--------------------|---|-------|-------|----------------------|-------|-------|----------------------|-------|-------|----------------------|-------|-------|----------------------|-------|-------|
| Temperature range | AE | BE | CE | AE | BE | CE | AE | BE | CE | AE | BE | CE | AE | BE | CE |
| | -40°C | -55°C | -86°C | -40°C | -55°C | -86°C | -40°C | -55°C | -86°C | -40°C | -55°C | -86°C | -40°C | -55°C | -86°C |
| External size | W1365 x D720 x H1101 | | | W1665 x D720 x H1101 | | | W1965 x D720 x H1101 | | | W2265 x D720 x H1101 | | | W2565 x D720 x H1101 | | |
| Internal size | W600 x D470 x H710 | | | W900 x D470 x H710 | | | W1200 x D470 x H710 | | | W1500 x D470 x H710 | | | W1800 x D470 x H710 | | |
| Inside measurement | 200L | | | 300L | | | 400L | | | 500L | | | 600L | | |
| Electric capacity | Capacity necessary for installment AC220V 1ph (50/60Hz) | | | | | | | | | | | | | | |
| Weight | 255Kg | | | 290Kg | | | 315Kg | | | 340Kg | | | 365Kg | | |

Options

Rack / Box / CO2 Backup system / 7"LCD touch screen and USB (3 years Data) - Basic -1channel (sensor 3channel, 6channel options) / U-system (remote alarm system and wireless control system – telecommunication subscription needed / Automatic voltage regulator / Cryogenic gloves / Battery backup (Control panel and alarm) / Recorder / VIP / ※The minimum temperature of ultralow temperature freezer(C-type) is supplied with -91°C performance if the customer requests.

Ultralow temperature freezer-DFU

Upright type Single door

Product features and specifications

- Realization of powerful and stable -86°C freezing performance with application of Two stage Auto cascade system
- The controller attached at the user's eye level
- Use of CFC-free special mixed refrigerants
- Application of aerodynamically designed pan wings equipped with huge capacity and high performance
- Hermetic compressor for LBP with excellent performance
- Use of 130mm high density urethane insulating materials to minimize cold air loss
- Non scaling stainless chamber and four stainless doors inside
- Strong structure frame and impact resistant powder coating finishing
- Round corner design for easy cleaning
- Prevention of condensation and the chill loss with three-layered wrinkle packing
- Filter cover for easy cleaning and washable filter
- Chest type cabinet structure
- Temperature for hard condition resistant test +32°C
- Temperature measuring sensor: platinum PT-100Ω(Class A 0.15 grade)
- Prevention of frost with Rim Heating system
- RS-232 Data port and Interface which is available for store, control and monitor temperature data via user's PC
- Equipped with the port for auto vacuum release



Control system

- Micro processor control system /Digital temperature indicator which adjust the temperature by 0.1°C
- Functions for setting double passwords to prevent access or control except administrators
- Data locking function

Alarm system and safety devices

- Alarm for filter cleaning period (Warm Condenser alarm)
- Alarm for high/low temperature of sight and auditory
- Equipped with pressure switch
- Equipped with overheating prevention switch
- Auto return system after blackout

Ultralow temperature freezer (Upright freezer)

| Model | DFU-128 | | | DFU-256 | | | DFU-374 | | | DFU-446 | | | DFU-558 | | | DFU-657 | | | DFU-740 | | | DFU-868 | | |
|--------------------|---|------|------|---------------------|------|------|---------------------|------|------|---------------------|------|------|----------------------|------|------|----------------------|------|------|-----------------------|------|------|-----------------------|------|------|
| Temperature range | AE | BE | CE | AE | BE | CE | AE | BE | CE | AE | BE | CE | AE | BE | CE | AE | BE | CE | AE | BE | CE | AE | BE | CE |
| | -40℃ | -55℃ | -86℃ | -40℃ | -55℃ | -86℃ | -40℃ | -55℃ | -86℃ | -40℃ | -55℃ | -86℃ | -40℃ | -55℃ | -86℃ | -40℃ | -55℃ | -86℃ | -40℃ | -55℃ | -86℃ | -40℃ | -55℃ | -86℃ |
| External size | W710 x D750 x H1390 | | | W710 x D750 x H1990 | | | W860 x D800 x H1951 | | | W860 x D900 x H1951 | | | W1010 x D900 x H1901 | | | W1010 x D990 x H1901 | | | W1010 x D1010 x H2001 | | | W1010 x D1151 x H2001 | | |
| Internal size | W450 x D475 x H600 | | | W450 x D475 x H1200 | | | W600 x D520 x H1200 | | | W600 x D620 x H1200 | | | W750 x D620 x H1200 | | | W750 x D730 x H1200 | | | W750 x D750 x H1300 | | | W750 x D891 x H1300 | | |
| Inside measurement | 128L | | | 256L | | | 347L | | | 446L | | | 558L | | | 657L | | | 731L | | | 868L | | |
| Electric capacity | Capacity necessary for installment AC220V 1ph (50/60Hz) | | | | | | | | | | | | | | | | | | | | | | | |
| Weight | 210kg | | | 300kg | | | 370kg | | | 380kg | | | 400kg | | | 420kg | | | 450kg | | | 500kg | | |

Ultralow temperature freezer (Power Plus model-Upright type)

- Upright type Power Plus model- the freezer appropriate for the surrounding temperature +40°C and over

| Model | DFU-128P | | | DFU-256P | | | DFU-374P | | | DFU-446P | | | DFU-558P | | | DFU-657P | | | DFU-740P | | | DFU-868P | | |
|--------------------|---|-------|-------|---------------------|-------|-------|---------------------|-------|-------|---------------------|-------|-------|----------------------|-------|-------|----------------------|-------|-------|-----------------------|-------|-------|-----------------------|-------|-------|
| Temperature range | AE | BE | CE | AE | BE | CE | AE | BE | CE | AE | BE | CE | AE | BE | CE | AE | BE | CE | AE | BE | CE | AE | BE | CE |
| | -40°C | -55°C | -86°C | -40°C | -55°C | -86°C | -40°C | -55°C | -86°C | -40°C | -55°C | -86°C | -40°C | -55°C | -86°C | -40°C | -55°C | -86°C | -40°C | -55°C | -86°C | -40°C | -55°C | -86°C |
| External size | W710 x D750 x H1390 | | | W710 x D750 x H1990 | | | W860 x D800 x H1951 | | | W860 x D900 x H1951 | | | W1010 x D900 x H1901 | | | W1010 x D990 x H1901 | | | W1010 x D1010 x H2001 | | | W1010 x D1151 x H2001 | | |
| Internal size | W450 x D475 x H600 | | | W450 x D475 x H1200 | | | W600 x D520 x H1200 | | | W600 x D620 x H1200 | | | W750 x D620 x H1200 | | | W750 x D730 x H1200 | | | W750 x D750 x H1300 | | | W750 x D891 x H1300 | | |
| Inside measurement | 128L | | | 256L | | | 347L | | | 446L | | | 558L | | | 657L | | | 731L | | | 868L | | |
| Electric capacity | Capacity necessary for installment AC220V 1ph (50/60Hz) | | | | | | | | | | | | | | | | | | | | | | | |
| Weight | 220kg | | | 310kg | | | 380kg | | | 390kg | | | 410kg | | | 430kg | | | 460kg | | | 510kg | | |

Options

Rack / Box / CO2 Backup system / 7"LCD touch screen and USB (3 years Data) - Basic -1channel (sensor 3channel, 6channel options) / U-system (remote alarm system and wireless control system – telecommunication subscription needed / Automatic voltage regulator / Cryogenic gloves / Battery backup (Control panel and alarm) / Recorder / VIP / ※The minimum temperature of ultralow temperature freezer(C-type) is supplied with -91°C performance if the customer requests.

Ultralow temperature freezer-DFU/CUS

Ultra Safe™(Twin Heart™) Freezer

Product features and specifications

- Realization of -90°C freezing performance with application of patented Single Stage Auto Cascade Two compressor systems
- When one compressor is broken down, another compressor maintains -80°C (with amb. +25°C/ +32°C)
- When exchanging or repairing the compressor, it does not need to cut off the power, and the samples do not have to be moved to another freezer, Also, welding is not necessary, so welding equipment is not needed. Vacuuming is not necessary, and the user can exchange the compressor within 30 minutes without trained engineers or special tools, so it can reduce the repair cost by 30~70% comparing to the previous product.
- The controller attached at the user's eye level
- Use of CFC-free special mixed refrigerants
- Application of aerodynamically designed pan wings equipped with huge capacity and high performance
- Hermetic compressor for LBP with excellent performance
- Use of 130mm high density urethane insulating materials to minimize cold air loss
- Non scaling stainless chamber and four stainless doors inside
- Strong structure frame and impact resistant powder coating finishing
- Round corner design for easy cleaning
- Prevention of condensation and the chill loss with three-layered wrinkle packing
- Filter cover for easy cleaning and washable filter
- Upright type canbinet structure
- Temperature for hard condition resistant test +32°C
- Temperature measuring sensor: platinum PT-100Ω(Class A 0.15 grade)
- Prevention of frost with Rim Heating system
- RS-232 Data port and Interface which is available for store, control and monitor temperature data via user's PC
- Equipped with the port for auto vacuum release



Control system

- Micro processor control system /Digital temperature indicator which adjust the temperature by 0.1°C
- Functions for setting double passwords to prevent access or control except administrators
- Data locking function

Alarm system and safety devices

- Alarm for filter cleaning period (Warm Condenser alarm)
- Alarm for high/low temperature of sight and auditory
- Equipped with pressure switch
- Equipped with overheating prevention switch
- Auto return system after blackout

Ultralow temperature freezer (Ultra safe (twin Heart)® freezer)

| Model (-90°C) | DFU-558CUS | DFU-657CUS | DFU-740CUS | DFC-200CUS | DFC-300CUS | DFC-400CUS | DFC-500CUS | DFC-600CUS |
|--------------------|---|----------------------|-----------------------|----------------------|----------------------|----------------------|----------------------|----------------------|
| External size | W1010 x D900 x H1901 | W1010 x D990 x H1901 | W1010 x D1010 x H2001 | W1365 x D720 x H1101 | W1665 x D720 x H1101 | W1965 x D720 x H1101 | W2265 x D720 x H1101 | W2565 x D720 x H1101 |
| Internal size | W750 x D620 x H1200 | W750 x D730 x H1200 | W750 x D750 x H1300 | W600 x D470 x H710 | W900 x D470 x H710 | W1200 x D470 x H710 | W1500 x D470 x H710 | W1800 x D470 x H710 |
| Inside measurement | 558L | 657L | 731L | 200L | 300L | 400L | 500L | 600L |
| Electric capacity | Capacity necessary for installment AC220V 1ph (50/60Hz) | | | | | | | |
| Weight | 410kg | 425kg | 455kg | 255Kg | 290Kg | 315Kg | 340Kg | 365Kg |

Options

1. Rack
2. Box
3. CO2 Backup system
4. 7"LCD touch screen and USB (3 years Data) - Basic -1channel (sensor 3channel, 6channel options)
5. U-system (remote alarm system and wireless control system – telecommunication subscription needed)
6. Automatic voltage regulator
7. Cryogenic gloves
8. Battery backup (Control panel and alarm)
9. Recorder
10. VIP

Ultralow temperature freezer-DFUD

Upright type double door with double controller

Product features and specifications

- Creative upright type double door freezer by OPERON which is independently controlled
- Independent temperature adjustment and space usage according to the characteristics of samples
- To keep samples which are used frequently in the upper chamber, and to keep samples which need long-term preservation in the lower chamber
- More stable temperature with triple silicon gasket of the upper and the lower chambers
- Four inner doors which are two in the upper chamber and two in the lower chamber and two outer doors in the upper and the lower chamber



Ultralow temperature freezer(Double door / Double controller freezer)

| Model | DFUD-558 | | | | | | DFUD-657 | | | | | | DFUD-740 | | | | | |
|----------------------------------|---|-------|-------|-------|-------|-------|-----------------------|-------|-------|-------|-------|-------|-----------------|-------|-------|-------|-------|-------|
| | AA | BB | CC | AA | BB | CC | AA | BB | CC | AA | BB | CC | AA | BB | CC | AA | BB | CC |
| Temperature in the upper chamber | -40°C | -55°C | -86°C | -40°C | -55°C | -86°C | -40°C | -55°C | -86°C | -40°C | -55°C | -86°C | -40°C | -55°C | -86°C | -40°C | -55°C | -86°C |
| Temperature in the lower chamber | -40°C | -55°C | -86°C | -40°C | 55°C | -86°C | -40°C | -55°C | -86°C | -40°C | -55°C | -86°C | -40°C | -55°C | -86°C | -40°C | -55°C | -86°C |
| External size | W1010 x D900 x H1990 | | | | | | W1010 x D1030 x H1950 | | | | | | Customised item | | | | | |
| Internal size | W750 x D620 x H550 | | | | | | W750 x D770 x H530 | | | | | | | | | | | |
| Inside measurement | 511L | | | | | | 612L | | | | | | | | | | | |
| Door (inside/outside) | 4/2 | | | | | | 4/2 | | | | | | 4/2 | | | | | |
| Electric capacity | Capacity necessary for installment AC220V 1ph (50/60Hz) | | | | | | | | | | | | | | | | | |
| Weight | 420kg | | | | | | 430kg | | | | | | | | | | | |

Upright type double door with single controller

It is an upright type double door freezer to minimize the cold air loss when taking out samples. As samples which are used frequently are kept in the upper chamber and samples which need a long-term preservation are kept in the lower chamber, it affects the temperature of other samples less.

Control system

- Micro processor control system /Digital temperature indicator which adjust the temperature by 0.1°C
- Functions for setting double passwords to prevent access or control except administrators
- Data locking function

Alarm system and safety devices

- Alarm for filter cleaning period (Warm Condenser alarm)
- Alarm for high/low temperature of sight and auditory
- Equipped with pressure switch
- Equipped with overheating prevention switch ■ Auto return system after blackout



Ultralow temperature freezer(Double door / Single controller freezer)

| Model | DFUD-374 | | | DFUD-446 | | | DFUD-558 | | | DFUD-657 | | | DFUD-740 | | |
|--------------------|---|------|------|---------------------|------|------|----------------------|------|------|----------------------|------|------|-----------------------|------|------|
| Temperature range | AE | BE | CE | AE | BE | CE | AE | BE | CE | AE | BE | CE | AE | BE | CE |
| | -40℃ | -55℃ | -86℃ | -40℃ | -55℃ | -86℃ | -40℃ | -55℃ | -86℃ | -40℃ | -55℃ | -86℃ | -40℃ | -55℃ | -86℃ |
| External size | W860 x D800 x H1951 | | | W860 x D900 x H1951 | | | W1010 x D900 x H1901 | | | W1010 x D990 x H1901 | | | W1010 x D1010 x H2001 | | |
| Internal size | W600 x D520 x H1200 | | | W600 x D620 x H1200 | | | W750 x D620 x H1200 | | | W750 x D730 x H1200 | | | W750 x D750 x H1300 | | |
| Inside measurement | 347L | | | 446L | | | 558L | | | 657L | | | 731L | | |
| Electric capacity | Capacity necessary for installment AC220V 1ph (50/60Hz) | | | | | | | | | | | | | | |
| Weight | 370kg | | | 380kg | | | 400kg | | | 415kg | | | 445kg | | |

Options

Rack / Box / CO2 Backup system / 7"LCD touch screen and USB (3 years Data) - Basic -1channel (sensor 3channel, 6channel options) / U-system (remote alarm system and wireless control system – telecommunication subscription needed / Automatic voltage regulator / Cryogenic gloves / Battery backup (Control panel and alarm) / Recorder / VIP / ※The minimum temperature of ultralow temperature freezer(C-type) is supplied with -91°C performance if the customer requests.

Ultralow temperature freezer

DFU-V / DFUD-V (shelf-cooling type)
DFU-P / DFUD-P (Power Plus model)

Upright type Validity(Shelf Cooling Type) – (DFU-V / DFUD-V)

Product features and specifications

- Realization of strong cooling power with installment of copper pipes in each shelf and accurate and even temperature profile inside



Ultralow temperature freezer (Shelf cooling type Validity freezer)-Single door

| Model | DFU-128V | | | DFU-256V | | | DFU-374V | | | DFU-446V | | | DFU-558V | | | DFU-657V | | | DFU-740V | | | DFU-868V | | |
|--------------------|---|------|------|---------------------|------|------|---------------------|------|------|---------------------|------|------|----------------------|------|------|----------------------|------|------|-----------------------|------|------|-----------------------|------|------|
| Temperature range | AE | BE | CE | AE | BE | CE | AE | BE | CE | AE | BE | CE | AE | BE | CE | AE | BE | CE | AE | BE | CE | AE | BE | CE |
| | -40℃ | -55℃ | -86℃ | -40℃ | -55℃ | -86℃ | -40℃ | -55℃ | -86℃ | -40℃ | -55℃ | -86℃ | -40℃ | -55℃ | -86℃ | -40℃ | -55℃ | -86℃ | -40℃ | -55℃ | -86℃ | -40℃ | -55℃ | -86℃ |
| External size | W710 x D750 x H1390 | | | W710 x D750 x H1990 | | | W860 x D800 x H1951 | | | W860 x D900 x H1951 | | | W1010 x D900 x H1901 | | | W1010 x D990 x H1901 | | | W1010 x D1010 x H2001 | | | W1010 x D1151 x H2001 | | |
| Internal size | W450 x D475 x H600 | | | W450 x D475 x H1200 | | | W600 x D520 x H1200 | | | W600 x D620 x H1200 | | | W750 x D620 x H1200 | | | W750 x D730 x H1200 | | | W750 x D750 x H1300 | | | W750 x D891 x H1300 | | |
| Inside measurement | 128L | | | 256L | | | 347L | | | 446L | | | 558L | | | 657L | | | 731L | | | 868L | | |
| Electric capacity | Capacity necessary for installment AC220V 1ph (50/60Hz) | | | | | | | | | | | | | | | | | | | | | | | |
| Weight | 210kg | | | 300kg | | | 370kg | | | 380kg | | | 400kg | | | 420kg | | | 450kg | | | 500kg | | |

Ultralow temperature freezer (Shelf cooling type Validity freezer)-Double door

| Model | DFUD-374V | | | DFUD-446V | | | DFUD-558V | | | DFUD-657V | | | DFUD-740V | | | DFUD-1080V | | | DFUD-1215V | | DFUD-1260V | |
|--------------------|---|-------|-------|---------------------|-------|-------|----------------------|-------|-------|----------------------|-------|-------|-----------------------|-------|-------|-----------------------|-------|-------|-----------------------|-------|-----------------------|-------|
| Temperature range | AE | BE | CE | AE | BE | CE | AE | BE | CE | AE | BE | CE | AE | BE | CE | AE | BE | CE | AE | BE | AE | BE |
| | -40°C | -55°C | -86°C | -40°C | -55°C | -86°C | -40°C | -55°C | -86°C | -40°C | -55°C | -86°C | -40°C | -55°C | -86°C | -40°C | -55°C | -86°C | -40°C | -55°C | -40°C | -55°C |
| External size | W860 D800 H1951 | | | W860 x D900 x H1951 | | | W1010 x D900 x H1901 | | | W1010 x D990 x H1901 | | | W1010 x D1010 x H2001 | | | W1460 x D1010 x H1990 | | | W1610 x D1010 x H1990 | | W2260 x D1010 x H1583 | |
| Internal size | W600 D520 H1200 | | | W600 x D620 x H1200 | | | W750 x D620 x H1200 | | | W750 x D730 x H1200 | | | W750 x D750 x H1300 | | | W1200 x D750 x H1200 | | | W1350 x D750 x H1200 | | W1400 x D750 x H1200 | |
| Inside measurement | 347L | | | 446L | | | 558L | | | 657L | | | 731L | | | 1080L | | | 1215L | | 1260L | |
| Electric capacity | Capacity necessary for installment AC220V 1ph (50/60Hz) | | | | | | | | | | | | | | | | | | | | | |
| Weight | 380kg | | | 390kg | | | 410kg | | | 425kg | | | 455kg | | | 500kg | | | 520kg | | 540kg | |

Upright-Double door type Power plus model

Product features and specifications

- The freezer appropriate for the surrounding temperature +40°C and over

Ultralow temperature freezer (Power Plus model-Double door Upright type)

| Model | DFUD-374P | | | DFUD-446P | | | DFUD-558P | | | DFUD-657P | | | DFUD-740P | | | DFUD-1080P | | | DFUD-1215P | | DFUD-1260P | |
|--------------------|---|-------|-------|---------------------|-------|-------|----------------------|-------|-------|----------------------|-------|-------|-----------------------|-------|-------|-----------------------|-------|-------|-----------------------|-------|-----------------------|-------|
| Temperature range | AE | BE | CE | AE | BE | CE | AE | BE | CE | AE | BE | CE | AE | BE | CE | AE | BE | CE | AE | BE | AE | BE |
| | -40°C | -55°C | -86°C | -40°C | -55°C | -86°C | -40°C | -55°C | -86°C | -40°C | -55°C | -86°C | -40°C | -55°C | -86°C | -40°C | -55°C | -86°C | -40°C | -55°C | -40°C | -55°C |
| External size | W860 x D800 x H1951 | | | W860 x D900 x H1951 | | | W1010 x D900 x H1901 | | | W1010 x D990 x H1901 | | | W1010 x D1010 x H2001 | | | W1460 x D1010 x H1990 | | | W1610 x D1010 x H1990 | | W2260 x D1010 x H1583 | |
| Internal size | W600 x D520 x H1200 | | | W600 x D620 x H1200 | | | W750 x D620 x H1200 | | | W750 x D730 x H1200 | | | W750 x D750 x H1300 | | | W1200 x D750 x H1200 | | | W1350 x D750 x H1200 | | W1400 x D750 x H1200 | |
| Inside measurement | 347L | | | 446L | | | 558L | | | 657L | | | 731L | | | 1080L | | | 1215L | | 1260L | |
| Electric capacity | Capacity necessary for installment AC220V 1ph (50/60Hz) | | | | | | | | | | | | | | | | | | | | | |
| Weight | 380kg | | | 390kg | | | 410kg | | | 425kg | | | 455kg | | | 500kg | | | 520kg | | 540kg | |

Control system

- Micro processor control system /Digital temperature indicator which adjust the temperature by 0.1°C
- Functions for setting double passwords to prevent access or control except administrators
- Data locking function

Alarm system and safety devices

- Alarm for filter cleaning period (Warm Condenser alarm)
- Alarm for high/low temperature of sight and auditory
- Equipped with pressure switch
- Equipped with overheating prevention switch
- Auto return system after blackout

Options

Rack / Box / CO2 Backup system / 7"LCD touch screen and USB (3 years Data) - Basic -1channel (sensor 3channel, 6channel options) / U-system (remote alarm system and wireless control system – telecommunication subscription needed / Automatic voltage regulator / Cryogenic gloves / Battery backup (Control panel and alarm) / Recorder / VIP / ※The minimum temperature of ultralow temperature freezer(C-type) is supplied with -91°C performance if the customer requests.

Accessory

Full Storage Capacity

Operon Chest freezer's Storage system(Rack & Box)

| Model | Eight stage | The number of box | Ten stage | The number of box | Twelve stage | The number of box |
|---------|-------------|-------------------|-----------|-------------------|--------------|-------------------|
| DFC-84 | 6EA | 48 | - | - | - | - |
| DFC-200 | 12EA | 96 | 12EA | 120 | 12EA | 144 |
| DFC-300 | 18EA | 144 | 18EA | 180 | 18EA | 216 |
| DFC-400 | 24EA | 192 | 24EA | 240 | 24EA | 288 |
| DFC-500 | 30EA | 240 | 30EA | 300 | 30EA | 360 |
| DFC-600 | 36EA | 288 | 36EA | 360 | 36EA | 432 |
| CFQ-150 | 8EA | 64 | 8EA | 80 | - | - |
| CFQ-152 | 8EA | 64 | 8EA | 80 | - | - |
| CFQ-156 | 8EA | 64 | 8EA | 80 | - | - |
| CFQ-232 | 15EA | 120 | 15EA | 150 | - | - |
| CFQ-300 | 21EA | 168 | 21EA | 210 | - | - |

Operon Upright freezer's Storage system(Rack & Box)

| Model | The number of shelf | Rack (stage X row) | 2"box+Cell divider | Micro tube (81 hole) | Micro tube (100 hole) |
|-------------|---------------------|-------------------------|--------------------|----------------------|-----------------------|
| DFU-128 | 3 | 9EA (3x3) | 81 | 6561 | 8100 |
| DFU-256 | 4 | 12EA (5x3) | 180 | 14580 | 18000 |
| DFU-374 | 4 | 16EA (5x3) | 240 | 19440 | 24000 |
| DFU-446 | 4 | 16EA (5x4) | 320 | 25920 | 32000 |
| DFU-558 | 4 | 20EA (5x4) | 400 | 32400 | 40000 |
| DFU-657 | 4 | 20EA (5x5) | 500 | 40500 | 50000 |
| DFU-740 | 4 | 15EA (5x5) + 5EA (6x5) | 525 | 42525 | 52500 |
| DFU-868 | 4 | 10EA (6x6) + 10EA (5x6) | 660 | 53460 | 66000 |
| DFUD-1080EV | 4 | 24EA (5x5) + 8EA (5x4) | 760 | 61560 | 76000 |



Upright rack



Sliding rack



Chest rack



Poly Carbonate Box (2")



Vacuum Insulation Panel (V.I.P)



U-System
(remote alarm, wireless control system)



Frost remover



Cryogenic gloves(Cryo gloves)



Circular 7 days chart recorder



Konics recorder(KRN-50)-2ch



Paperless recorder(2ch~12ch)




Yokogawa recorder SR-10006(1~6ch)

Description

Controller 1 (LCD controller: option)

Dual cooling system which does not need an oil separator (Auto cascade system)



LED controller (standard)

- Data Lock function: Set-point security password
- The temperature is set by 0.1°C, and setting temperature and low/high temperature alarm can be entered. Control circuit chooses safety circuit when operated again after blackout to protect the device in case of auto power return.
- When choosing remote alarm system, alarm messages about abnormal temperature and blackout are sent to maximum four users.
- RS-232 PC communications program is available for monitoring and setting temperature such as accumulated time/temperature recording/graph presentation.
- Text description and LED flickering while sending auditory alarm are presented about abnormal temperature and abnormal condensation.



| START/STOP | Start/stop |
|--------------|---|
| MODE | Save data or Move to the next stage |
| | High temperature alarm/Low temperature alarm/Remote alarm |
| | LN2.CO2 back up alarm |
| | Data lock, set point security |
| | Setting password etc |
| SHIFT | Move the digit |
| INC | Increase of numbers |
| ALARM ON/OFF | Alarm ON/OFF |

New controller

- 7"Color TFT-LCD, Touch Screen Panel
- Setting the temperature by 0.1°C and entering setting temperature and low/high temperature alarm, set-point security and password security
- Control section(Touch screen panel) includes data acquisition function, so additional recorder devices are not necessary.
- Data acquisition and backup of collected data using USB are easy, and save them as an excel file.
- Validation Support System presents maximum six temperature and graphs (including control sensor) at the same time and is available for setting save period from per a minute to per an hour recording maximum 3 years. Also, real-time monitoring is available performing the function of validator (OQ, PQ, Paperless Recorder, etc). (Deluxe Pack)
- (Option for additional sensors), 1-(Basic Pack control 1), 2-(Ice Crystal Pack control1, ice crystal temperature of samples monitoring 1, Abnormal condensation 1), 3-(Deluxe Pack control1, shelf3, abnormal condensation1, outdoor temperature1)
- Real-time graph can be checked on graph monitoring mode, so control state of products can be checked at a glance.
- Energy saving function is set to realize energy saving
- Visual alarm screen is made, so various alarm and warning messages are presented on Touch Screen Panel in the form of graphics or texts. Information on abnormal high temperature, low temperature or blackout is sent to the user.
- If UPS is chosen, it supplies the power to control section and sends information on blackout and various alarms to the user. Saving data during blackout is available. Automatic return to already saved set value in case of power return after blackout. (When choosing U-System, alarm messages are sent to the mobile)

LCD(7") touch screen controller (option)



(It may be different from the actual image)


Controller & Alarm system – LCD touch screen control panel (option)

- Ultralow temperature freezer communicates with PC using U-system through RS-485 in the remote area which make monitoring, control and alarm available. One system controls maximum twelve ultralow temperature freezers.
- CDMA of U-system makes the user control twelve ultralow temperature freezers, monitor the current state of the device, receive alarm message(text message) and backup accumulated alarms with the mobile by wireless realizing safer sample management.

Code classification

Temperature Code

| CODE | TEMPERATURE |
|------|---------------|
| A | -40°C |
| B | -55°C |
| C | -86 °C |
| D | -120°C |
| CRYO | -150°C~-203°C |



Model Code

| Code | Model | REMARKS |
|----------|--|---|
| DFU | Deep Freezer Upright type | DFU-128/256/374/446/558/657/740L |
| DFU-V | Deep Freezer Upright type with Validation | Plumbing each shelves |
| DFUD-V | Deep Freezer Upright type Double door with Validation | DFU-374/446/558/657/740/1080/1260L |
| DFUD-E | Deep Freezer Upright type Double door | Single controller(same as DFU) |
| DFUD | Deep Freezer Upright type Top chamber:-86°C/-55°C/-40°C Bottom chamber:-86°C/-55°C/-40°C | Double controller DFUD-558CC/558AC/558BB/558BC/558AA DFUD-657CC/657AC/657BB/657BC/657AA |
| DFUP | Deep Freezer Upright Power plus up | Frequently door open or poor environment |
| DFC | Deep Freezer Chest type | DFC-84/200/300/400/500/600L |
| CFQ | Cryo Freezer Conqueror | CFQ-150/-152/-156/-300/-232 |
| FDB | Freeze Dryer Bench top type | FDB-5503 |
| FDU | Freeze Dryer Upright type | FDU-7003/7006/7012/7024/8603/8606/8612/8624 |
| FDG | Freeze Dryer Glass type | FDG-120/105/90 |
| FDS/FDTS | Freeze Dryer Stoppering type | FDS-12012/12006/8612/8606/7003/5503 , FDTS(For Production) |
| FDT | Freeze Dryer Bulk Tray type | FDT-8611/5511/8650/8632/8612/8606/12050/12032/12012/12006 |
| FDCF | Freeze Dryer Chemical Free | FDCF-12012/12006/12003 |
| SVQ | Speed Vacuum concentrator | SVQ-120/-95, MSVQ-20 |
| IMC | Immersion Cooler | IMC-120/-86 |
| SF | Shell Freezer-Cryo bath | SF-40/CB-80 |
| WVCP | Water Vapor Cryo Pump | WVCP-1500/1650 |
| CT | Cold Trap (Glass type/Coil type/Probe type) | CT-120/-90/-70/-55 , CTL/CTM/CTI |
| CH | Industrial Ultra Low Temp.of Chiller | CH-40/-55/-86/-120 |
| CTC | Cryogenic Test Chamber | CTC-165 /MTC |
| GR | Gas Recycler | GR-120/86/40/20 |
| CnH | Circulation Bath(cooling & heating) | CnH-2505/2510/2520/2530/2540/2550/4005/4010/4020/4030/4040/4050 |
| HC | Circulation Bath (heating) | HC-20005/20010/20020/20030/20040/20050 |
| IC | Industrial chiller | IC-4020/6020/9020/4035/6035/9035 |
| FDUT | Freeze Dryer (Combi) | FDUT-12012/12006/12003/8612/7012/8606/7006/8603/7003 |
| MPS | Multi Purpose System | MPS-55 |
| CUS | Ultra Safe (twin heart) Freezers | DFU-CUS/DFUD-CUS/DFC-CUS |
| ATC/VTC | Altitude-Temperature-Vacuum test chamber | ATC-100/200/300 , VTC-100/200/300 |


Controller Code

| Code | Type | Model | Remarks |
|------|-------------------|--|------------------|
| C/E | Compact/Economic | All kinds of freezer DFC-84/IMC/CT/SF-40 | For freezer |
| E | Economic Standard | Basic control panel for freeze dryers | For Freeze Dryer |
| P | Programmable | FDT-additional drying program | |
| N | Premium Network | All Models(option), Touch screen with Memory card & Tele-communication | |

Application Code

| Classification | Model | Remarks |
|--|----------------------------------|------------|
| Freezer for sample preservation | DFC/DFU/DFUD/DFUV/DFUP/CUS | |
| Epoxy freezer / Ice pack freezer | DFUD-1080EV/1260EVtype | |
| Plasma Freezer | DFU-446/CFQ-156 | |
| STEM CELL Freezer | CFQ | |
| Compact Freezer for personal use | DFC-84 | |
| Freeze Dryer for laboratory | FDB/FDU/FDT/FDCF/FDG/FDS/MPS | |
| Freeze Dryer for food and medicines (Cooled & Heated shelf) | FDT, FDUT | |
| Freeze Dryer for collect of organic solvents like ethanol | FDG/FDCF | |
| Cryo pump for semiconductor/Sputtering/Coater | WVCP-1500/WVCP-1650 | Order made |
| Concentrator for analysis of DNA/RNA/PROTEINS | SVQ, MSVQ-20 | |
| Cooler for a partial cooling | IMC | |
| Ultra-low chiller for industrial use | CH-40/-55/-86/-120 | Order made |
| Quick Freezer | LT-40 | Order made |
| Cryogenic test chamber(Strength) | CTC-165/150/130 , MTC-150/130 | Order made |
| Circulation Bath for a reaction unit or distillation unit, chiller | CnH/HC/CH/IC | Order made |
| Collecting Vapor, Gas, Organic Solvent | CT-50/70/90/120, CTL/CTM/CTI | |
| Altitude-Temperature-Vacuum Test chamber | ATC-100/200/300, VTC-100/200/300 | Order made |

Ultra-Low Temperature Range



Vacuum Relief Port-Built in
Rim Heating Systems-Built in

| Type | Data Setting | | | | | | | | Alarm | | | | | | MEMORY CARD |
|--------|-----------------|-------------|-------------|-----------------|---------------------|---------|-----------|-----------------|-----------------|---------------|----------------|--------------|------------|--------|-------------|
| | MICRO PROCESSOR | TEMPERATURE | ALARM TEMP. | ACCESS PASSWORD | DATA RESET PASSWORD | DATALOC | PRINT SET | C.02 BACKUP SET | BUZER& LED LAMP | POWER FAILURE | WARM CONDENSER | CLEAN FILTER | ROOM TEMP. | R.T.M. | |
| E-Type | Yes | Yes | Yes | Yes | Yes | Yes | Yes | Yes | Yes | OPTION | OPTION | OPTION | OPTION | OPTION | - |
| N-Type | Yes | Yes | Yes | Yes | Yes | Yes | Yes | Yes | Yes | Yes | Yes | Yes | Yes | Yes | Yes |