

SiriusWeb Pharmacy

Monitor fridge temperatures from the cloud



Cold chain compliance is becoming more subjected to external regulation. Sirius Web is a cloud based solution allows on-line access to review your fridge and medication temperatures from any web or iPhone connection.

- Measures and monitors your fridge storage conditions
- Automatically records and archives all your temperature data
- On-line access for review available from any PC connected to the internet
- Automatic alarm alerts sent via email, sms/text message or phone
- Easy to install and use
- User friendly web interface
- iPhone App

Easy solution for MOH compliance
No manual data download
24/7 monitoring
Peace of mind protection of
expensive medicines and vaccines



Available on the
App Store

Check your
temperature data
on your iPhone

SYSTEM ARCHITECTURE

Step 1 : Install the data logger on the fridge to be monitored.



The SPY RF recorder measures and records the temperature data and wirelessly transmits to the Sirius Box receiver.

Step 2 : Position the SiriusBox on wall and simply plug into a power outlet and internet.




The Sirius Box then transfers measured data daily through internet the Sirius Web server.

Step 3 : Review data and share information through any PC or web device on Sirius Web site.



In case of alarm, the recorder sends an alert message to the Sirius Web server via the SiriusBox and pre-set alarm notifications such as email, text message or voice message can be launched.

- Unlimited access to Sirius Web
- Automatic backup of your data for up to 10 years (depending on package)
- Automatically emailed reports
- Local Alarm Buzzer
- Fridge Power Failure Alarm
- 2 year Plan
- Monitor up to 20 fridges and medication rooms
- Calibration service option



Monitoring is our priority